

# Use LabTemp® 20 To Monitor The Temperature Profile of Your Samples



Remove refrigerated sample from the refrigerator




**20 RED/17 BLUE**  
Sample is ready to use



**20 is GREEN**  
and sample is ready  
to be returned  
to the refrigerator



SHOWN ACTUAL SIZE  
(at refrigerated temperature)



COLOR INTERPRETATION

- TAN/RED:**  
Temperature is 1°C lower than the digit that appears
- GREEN:**  
Temperature is indicated by the digit that appears
- BLUE:**  
Temperature is 1°C higher than the digit that appears

