

LabTemp® 20 Indicator provides instant temperature monitoring at a glance.

Most of the leading manufacturers of biological compounds/chemicals have warnings on their labels to keep their material stored at specific temperatures (generally 2°- 8°C or be kept frozen). As most technicians well know, not only is it important to maintain these storage temperatures but also once out of storage to allow these materials **to equilibrate prior** to use, such as pipetting. LabTemp® 20 reversible temperature indicators can be affixed directly on the container before or after removal from storage and used as a thermometer to let the technician know when the **material has reached its equilibration temperature** and to monitor the materials while out of storage so they do not exceed specified temperatures.



2	BSI
5	20
8	17
11	14

Color Interpretation

TAN 5	GREEN 5	BLUE 5
5	5	5
means 4°C	means 5°C	means 6°C

LabTemp® 20 indicators are general purpose temperature indicators which accurately monitor temperatures between 1-21°C. LabTemp® 20 indicators should be used on all temperature sensitive materials used in laboratories **especially when these materials are removed from or need to be returned to a refrigerated environment.**

