



BIOSYNERGY, INC.

1940 E. Devon Avenue, Elk Grove Village, Illinois 60007
(847) 956-0471 Toll Free (800) 255-5274 Fax (847) 956-6050
e-mail: biosynergy@att.net website: www.biosynergyinc.com

LABTEMP® 20 SURFACE TEMPERATURE INDICATORS

LabTemp® 20 indicators are general purpose temperature indicators which accurately monitor temperatures between 1-21°C. LabTemp® 20 indicators should be used on all temperature sensitive materials in laboratories, **especially when these materials are removed from or need to be returned to a refrigerated environment.**

Most of the leading manufacturers of biological compounds/chemicals have warnings on their labels to keep their material stored at specific temperatures (generally 2°C – 8°C or be kept frozen). Not only is it important to maintain these storage temperatures, but also once out of storage it is essential to allow these materials **to equilibrate prior** to use. LabTemp® 20 reversible temperature indicators can be affixed directly on the container before or after removal from storage to be used as a thermometer to let the technician know when the **material has reached its equilibration temperature** and to monitor the materials while out of storage so they do not exceed specified temperatures.

LabTemp® 20 liquid crystal disposable indicators are designed for use in clinical chemistry, microbiology, and any laboratory area when it is critical that compounds and organisms are kept at specific temperatures. They are inexpensive disposable indicators and may be affixed to all types of plastic, glass, metal, and wood labware as well as curved, round and flat surfaces for up to 18 months.

LabTemp® 20 can be used on items stored in refrigerators, freezers, temperature controlled centrifuges, low temperature incubators, and even water baths for short term immersion in water (1-2 hours). This indicator provides instantaneous and continuous monitoring of your temperature sensitive materials at a glance.

The necessity for temperature monitoring of sensitive materials is a must in today's laboratory! Biosynergy's LabTemp® 20 indicators can be an effective, easy to use way of monitoring materials to meet storage requirements and **to make sure the materials have met their proper equilibration temperature for accuracy in testing.**

For a free sample or additional information, please call toll-free, 800-255-5274, or visit our website at www.biosynergyinc.com.