


**CONTROL ENZYME LOSS IN BLOOD SAMPLES**

Use **LabTemp®20** To Control Temperature Base

2	BSI
5	20
8	17
11	14


LabTemp® 20 Indicator

2.0°C  
(2 is GREEN)



LabTemp® 20

**EFFECT OF TEMPERATURE**



*Typical progress curves of an enzyme reaction at various temperatures, showing dependence of apparent optimum temperature on time.*

**COLOR INTERPRETATION**

- TAN/RED:**  
Temperature is 1°C lower than the digit that appears
- GREEN:**  
Temperature is indicated by the digit that appears
- BLUE:**  
Temperature is 1°C higher than the digit that appears

