

Use LabTemp® 20 To Monitor The Temperature Profile of Your Samples



Remove refrigerated sample from the refrigerator



20 RED/17 BLUE
Sample is ready to use



20 is GREEN
and sample is ready
to be returned
to the refrigerator



SHOWN ACTUAL SIZE
(at refrigerated
temperature)



COLOR INTERPRETATION

TAN/RED:
Temperature is 1°C lower than
the digit that appears

GREEN:
Temperature is indicated by
the digit that appears

BLUE:
Temperature is 1°C higher than
the digit that appears

